



# Amigo Classic FD

*Innovative Engineering  
Unlike Any Other!*

CLASSIC FD



Shown with  
SmartSensor  
Safety Option  
(Sold separately)

**NEW!**  
Rear Cover with  
Lifting Handle

The **Amigo Classic FD** • Item No. 113001  
**MSRP \$1994** with Batteries

The Amigo Classic FD is today's success version  
of the Front Drive first introduced in 1968.

- Functionality
- Durability
- Versatility
- Maneuverability
- Transportability

# FEATURES & BENEFITS


- Automatic Electromagnetic Brake with Manual Handbrake option for optimal safety
- 24-Volt Electronics offer highly efficient performance
- Infinitely Adjustable Handle assures steering-position comfort
- Narrow and compact 20" (50.8 cm) Rear Wheel Width
- 33" (83.8cm) Tight Turning Radius assures easy maneuverability
- Non Take-Apart Steel Platform defines rugged durability
- 250 Pounds (113.4 kg) robust weight capacity offers generous comfort
- Dynamic and Regenerative Braking
- Flip-up Adjustable, Removable Arms
- Fully Cushioned & Rotational Seat of 360° to position exact setting of choice
- Power Seat Lift Option maximizes user functionality

*The **Classic FD** is the perfect solution for your mobility needs.*



- NEW!**
- SmartSensor Safety Option prevents pedestrian impact
  - Infinitely Adjustable Seat Lock for select position
  - HI/LO Speed Controls
  - Oval Handle design protects control levers with easy grip perfect for limited hand usage

- NEW!**
- Color Selection: Midnight Black
  - Battery /Travel Distance:  
(2) 12AH / 18 miles  
(29 km)
  - Speed Range: 0-3.7mph  
(0-6 kph)
  - Flat-Free Front & Rear Tires
  - Programmable Controller
  - 2AMP-24V Offboard Charger w/AC Cord
  - Amigo limited lifetime/  
two-year warranty



**NEW!** Rear cover with lifting handle standard

YOUR LOCAL REPRESENTATIVE



Product pricing, options and accessories are subject to change without notice.