



Are you ready to go?

Amigo Mobility International myamigo.com • 1-989-777-0910



The Amigo HD is designed to provide a mobility solution for larger-framed individuals. With a solid steel frame and extended platform, the HD offers a unique combination of comfort and durability.



Large integrated seat pocket

FEATURES

Model	650000
Travel surfaces	Indoor/outdoor surfaces
Drive system	Heavy duty Unipac $^{\circ}$ &
	Vari-speed drive system
Platform	Solid steel, extended
Weight capacity	500 lbs (227 kg)
Speed	0-5 mph (0-8 kph)
Batteries	(2) 33/35 AH
	airline-safe, sealed
Battery capacity	30 miles (48.3 kilometers)
Charger	Onboard
Seat	22" (55.9 cm) vinyl folding seat
Front tire	8" (22.8 cm) foam filled
Rear tires	10.5" (26.6) foam filled
Turning radius	34" (96.5 cm)
Underbody clear	ance 3.4" (8.9 cm)
Length	45.5" (121.3 cm)
Rear wheel width 23" (59.7 c	
Height	37.5" (100.3 cm)
Total weight (with batteries) 170 lbs (84.4 kg)	

Extended seat lock with bend



Amigo Mobility International, Inc. Improving Lives Through Mobility® since 1968





*Product pricing, options and accessories are subject to change without notice.