

merits

ATLANTIS

K0826/K0827

HEAVY DUTY



FEATURES

- Weight Capacity 600 lbs.
- Heavy-duty power wheelchair with a unique and prominent design
- Four seat posts provide a more stable and durable structure for the frame
- Available with a pan or captain seat in 24 or 26-inch width
- Extra-wide seat and extra-high back provide great comfort for bariatric users
- Super high power in-line motors help riders negotiate 10° uphill ramps easily

SPECIFICATIONS

Weight Capacity	600 lbs
Max Speed	5 mph
Range	32 miles
Turning Radius	31.5"
Ground Clearance	3"
Overall Length	44"
Overall Width	32"
Batteries	Group 24 (sold separately)
Drive Wheels	12" Pneumatic tires
Casters	9" Front filled
Seat Width	22"-26"
Seat Depth	20"-22"
Back Height	20"
Controller	LinX
Warranty	12 months : Motor & Electronics 5 Year/ Limited Lifetime Frame

OPTIONS

Swing-away joystick
USB charger
Elevating legrests
Oxygen tank holder
Cane/Crutch holder
Walker holder
Foldable rear basket
Canopy
Height adjustable arms
Weather cover

MODELS

Atlantis w Captain Seat	P710-2VHD C
Atlantis w Pan Seat	P710-2VHD S



Merits Health Products



Meritshealthusa



Merits USA

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

239-239-0579 | sales@meritsusa.com

merits