

merits

DUALER

FULL SIZE



FEATURES

- Unit changes from rear to front-wheel drive with the swivel of a seat
- 5" power seat lift
- 16", 18" or 20" ergonomic semi-reclining seats; up to 135°
- 10" knobby drive tires for superior traction
- Dual footplates for comfort and convenience
- Red and blue color options

SPECIFICATIONS

Weight Capacity	300 lbs
Max Speed	5 mph
Per-Charge Range	18 miles
Batteries	U1-34Ah x 2
Ground Clearance	2"
Turning Radius	21"
Overall Length	38"
Overall Width	24"
Drive Wheels	10" foam filled tires
Casters	8" foam filled
Seat-to-Floor Height	21.5"-26"
Seat Width	16", 18" or 20"
Seat Depth	16"-18"
Back Height	18" or 20"
Controller	PG VR2
Power Seat Lift	5" - single post
Warranty	18 Months : Motor & Electronics 5 Year/ Limited Lifetime Frame 1 Year/ Batteries

OPTIONS

- Swing-away joystick
- USB charger
- Oxygen tank holder
- Foldable rear basket
- Cane & crutch holder
- Walker holder
- Weather cover
- Elevating legrests

MODELS

Dualer - FWD/RWD, Captain	P312-2C
Dualer - FWD/RWD, Pan	P312-2S



-  Merits Health Products
-  Meritshealthusa
-  Merits USA

Merits Health Products Inc.
 4245 Evans Avenue, Fort Myers FL, 33901
 Phone: 1-800-963-7487 Fax: 239-772-3252
 sales@meritsusa.com

