



# FR168-4S

**POPULAR  
PRACTICAL  
AFFORDABLE**



## SPECIFICATION

LENGTH	52.13 in (132.4 cm)
WIDTH	24.13 in (61.3 cm) at wheels
WEIGHT (SCOOTER)	146.74 lbs (66.7 kg)
WEIGHT (SCOOTER+BATTERIES)	214.94 lbs (97.7 kg)
WHEELS/TIRES	PU solid 10 in (25.4 cm)
TURNING RADIUS	62.6 in (159 cm)
GROUND CLEARANCE	2.6 in (6.6 cm)
SEAT	Foldable standard seating
MAX SPEED	4.97 miles/hr (8 kph)
BATTERY	2 x 50 Ah
WEIGHT CAPACITY	350 lbs (158.8 kg)
DRIVE SYSTEM	Rear-wheel drive, sealed transaxle, 24V DC motor
SUSPENSION	Front and rear shock absorber

The information contained herein is correct at the time of publication. Factors such as battery charge, battery condition, user weight, tire pressure, and terrain type will cause variation to maximum range and speed. Due to continual product improvement and manufacturing tolerances; length, width, turning radius and ground clearance can be subject to variance of +/- 3%. We reserve the right to alter specifications without prior notice.

## Shorter • Lighter • Steadier

A rear-view mirror, front and rear lighting, and turning indicators for safer travel in dimmed conditions.



Toll Free: (888) 957-1888 | Se Habla Español | 我們說中文  
Email: info@freeriderusa.com | Website: www.freeriderusa.com  
8696 Utica Ave, Rancho Cucamonga, CA 91730