

Geo Cruiser, portable mobility for everyday use,
and the explorer in you.

pathway
mobility



Features and Benefits

- + Lightweight and foldable design
- + Easy to transport - travel friendly
- + Safe and reliable - FDA cleared
- + Multi-terrain
- + Fits through standard home/business doorways
- + Armrests can be Raised for Easy Side Access
- + Removable seat covers for easy cleaning
- + Anodized anti-corrosion finish
- + Requires very little service and maintenance
- + Convenient under-seat basket



A lighter and easier way to get around.

The Geo Cruiser™ is a lightweight and foldable power chair that makes getting around easier and lighter.

It's portable and easy-fold design allows for convenient storage and easy transport. In fact, it fits in just about any vehicle, is approved for air travel, and is highly recommended by cruise lines.

Whether your customers are going shopping, to museums, amusement parks, family gatherings,

traveling or just wanting to gain more independence, they will experience more with the Geo Cruiser.

The Geo Cruiser is designed around an aircraft grade aluminum alloy frame, and is powered by two ultra-light lithium ion batteries and a set of twin brushless motors.

Navigating the Geo is a breeze using a quick-release multi-function joystick controller that can be easily programmed to meet a user's specific needs.

To personalize comfort, we included adjustable leg-guards, footrest, and a four position back-rest. The result, a machine that performs flawlessly with total comfort in mind, making it an ideal solution for everyday use as well as travel.



General Specifications

Geo DX

GEO LX

GEO ELITE EX

Product Weight: (with battery)	56 lbs.	59 lbs.	65 lbs.
Product Weight: (without battery)	50 lbs.	53 lbs.	59 lbs.
Maximum Weight Capacity:	270 lbs.	330 lbs.	350 lbs.
Unfolded Size: (L x W x H)	36" X 24.5" X 35"	36" X 26.5" X 35"	37" X 27.25" X 36.5"
Folded Size: (L x W x H)	13" X 24.5" X 28"	13" X 26.5" X 28.5"	13" X 27.25" X 30"
Travel Distance:	Up to 15 Miles	Up to 15 Miles	Up to 15 Miles
Maximum Speed: (adjustable)	Up to 4 mph	Up to 4 mph	Up to 5 mph
Seat Size: (L x W)	17" X 15"	17" X 17"	17" X 17"
Seat Height (from ground):	19.5"	20"	21"
Arm Rest Height (from ground):	28"	28"	28.5"
Turn Radius:	31.5"	34.5"	34.5"
Incline/Slope Capability:	12 Degrees	12 Degrees	12 Degrees
Motor:	(2) 180W Brushless Motors	(2) 180W Brushless Motors	(2) 250W Brushless Motors
Battery:	(2) Lithium Ion/24V 16Ah	(2) Lithium Ion/24V 16Ah	(2) Lithium Ion/24V 16Ah
Battery Weight:	6 lbs. (3 lbs. each)	6 lbs. (3 lbs. each)	6 lbs. (3 lbs. each)