






Standard Fabrics:




Calypso (TCP)



Camel (TCM)



Hazelnut (THA)



Shiraz (TSZ)



Sterling (TST)

This fabric stands out for its comfort, smooth texture, and rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

* An upcharge will apply for all crossover fabrics.

EZ Sleeper Slim PR763-MED



Comfort Zone™
by GOLDEN

Experience comfort like never before with the EZ Sleeper Slim by Golden. This model offers a modern style track arm design, while the plush back, chaise pad and supportive pocketed coil springs provide hours of relaxing comfort. The EZ Sleeper Slim comes standard with our 2-Zone Comfort Zone positioning system, which features our patented MaxiComfort® positioning operated by our programmable AutoDrive™ hand control.

Standard Features:

- Patented MaxiComfort® Positioning
- Programmable AutoDrive™ Hand Control
- Extra Pocket



MAXICOMFORT® Series PR763

Model Number	PR763-MED
Size	Medium
Class II Medical Device	Yes
Weight Capacity	375 LB
Positioning	MaxiComfort
Comfort Zone(s)	2
Twilight Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Bracket
Head Rest Cover	No
Arm Covers	No
Overall Width	33"
Overall Height	45"
Width Between Arms	23"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	30"
Seat to Top of Arm	8"
Distance Required from Wall Reclined	23"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Outside Arm
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
Heat Wave Technology Option	No
Articulating Headrest Option	No
Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Yes
USB Charging Port	No
Battery Backup	Yes
SmartTek® Diagnostic System	Standard
Place of Origin	USA



Shown in Brisa® Cream*

SmartTek® Diagnostic System STANDARD on all Golden models.
All Golden lift chairs are Class II Medical Devices.



Comfort Zone™
by GOLDEN

Comfort Zones 1 & 2 work together to create individualized wellness positioning and comfort. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

SHT-PR735-071723