

Standard Fabrics:



Acorn



Graphite



Maple

This faux leather is a rich and modern alternative to natural leather that is extremely comfortable while offering the rugged durability of leather. Made from 100% polyester, the fabric is easy to care for ensuring functionality as well as elegance. This long-lasting fabric will complement a variety of decors while standing the test of time under heavy use.

Crossover fabrics are unavailable.



Rhea PR442





Take a vacation from your busy day in The Rhea. It offers an overstuffed seam back design, plush chaise pad and an open arm design that holds numerous standard features that include two cup holders, two storage compartments, and a USB port. Positioning options are endless with our 4-Zone Comfort Zone™ Positioning System, which includes our patented MaxiComfort® Positioning, adjustable headrest and lumbar all operated by our programmable AutoDrive™ hand control. For extra wellness and comfort, our HeatWave™ Technology comes standard providing therapeutic warmth starting at your shoulders and continuing down to your feet. When you open the right arm, you will find a cup holder and a storage compartment that also holds the AutoDrive hand control. In the left arm, you will find an additional cup holder, a USB charging port, and a second storage compartment that holds the HeatWave hand control. As a bonus, you can use either cup holder without losing the comfort of the armrests thanks to the split design!

Standard Features:

- Patented MaxiComfort® Positioning
- Adjustable Headrest
- Adjustable Lumbar
- Programmable
 AutoDrive™ Hand Control
- HeatWave™ Technology
- USB Charging Port
- Built-in Cup Holders
- SmartTek® Diagnostic System
- Extra Pocket





DeLUNA Series PR442

| Model Number | PR442-MED | PR442-LAR |
|--|-------------|-------------------------|
| Size | Medium | Large |
| Class II Medical Device | YES | YES |
| Weight Capacity | 375 LBS | 375 LBS |
| Positioning | MaxiComfort | MaxiComfort |
| Comfort Zone(s) | 4 | 4 |
| Twilight Technology | NO | NO |
| Chaise Seat | YES | YES |
| Pocketed Coil Seat | YES | YES |
| Back Type | Seam | Seam |
| Removable Back Type | Bracket | Bracket |
| Head Rest Cover | NO | NO |
| Arm Rest Covers | NO | NO |
| Overall Width | 37" | 39" |
| Overall Height | 45" | 50" |
| Width Between Arms | 22" | 24" |
| Floor to Top of Seat | 20" | 21" |
| Seat Depth | 21" | 22" |
| Seat to Top of Back | 29" | 33" |
| Seat to Top of Arm | 7" | 8" |
| Distance Required from Wall Reclined | 25" | 30" |
| Extension Reclined from Back of Seat to End of Foot Rest | 43" | 44" |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | YES | YES |
| Warranty on Electrical Parts | 3 YR | 3 YR |
| Pro Rated Years 4 Thru 7 | YES | YES |
| Warranty on Mechanical Labor | 3 YR | 3 YR |
| Warranty on Electrical Labor | 1 YR | 1 YR |
| Hand Control Routing | Right Arm | Right Arm |
| | | Compartment |
| Extra Pocket Option | Standard | Standard |
| HeatWave Technology | Standard | Standard |
| Heat & Vibrating Massage Option | NO | NO |
| Articulating Headrest Option | Standard | Standard |
| Adjustable Lumbar Option | Standard | Standard |
| AutoFlex Option | NO | NO |
| Foot Rest Extension Option | NO | NO |
| Straight Lift Option | YES | YES |
| USB Charging Port | Left Arm | Left Arm |
| Battery Backup | Lithium-Ion | Compartment Lithium-Ion |
| Place of Origin | China | China |
| riace of Origin | Cilila | Cillia |

 ${\it SmartTek\ Diagnostic\ System\ STANDARD\ on\ all\ Golden\ models}.$

 $Designed \ and \ engineered \ in \ Old \ Forge, \ PA. \ Assembled \ in \ China \ using \ Golden \ components.$





STANDARD!

Far infrared heat penetrates skin and muscles to provide a soothing heat therapy experience. Adjustable Hi-Lo settings and timers let you control the comfort!

The Rhea comes standard with the following comfort and convenience features:



Right Arm: Cup Holder & Storage Compartment



Left Arm: USB Charging Port, Cup Holder & Storage Compartment



Adjustable Headrest & Lumbar



HeatWave™ Technology





Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.





UPPER BODY







TWILIGHT TILT

WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374 Fax: 800-628-5165 www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

ake changes to the product of illerature at any time.



