

Patriot[™]
4 wheel
GR575

Body Colors Red or

If you're looking for the ultimate outdoor driving experience, look no further than the Golden Patriot™. The Golden Patriot's rugged outdoor design features an off road package, four-wheel independent suspension, a full lighting package and a full sized captain's seat. With a 400 lb. weight capacity, 13" pneumatic tires, and a battery range of 24 miles, the Golden Patriot will give you the ride of your dreams.





Patriot 4 wheel **Specifications**

•	
Model	GR575
Medicare Code	Not Coded
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	8 mph
Operating Range ¹	24 miles
Ground Clearance	6.75"
Turning Radius	76.5"
Type Batteries	2-75 AH
Colors	Red, Orange
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Mechanical Hand Brake	Yes
Charger	Off Board
Controller	140 Amp
Dimensions:	
Length	57"
Width	28"
Height (ground to top	
of back of seat)	47" - 49"
Height (ground to top	
of back of headrest)	54.5" - 59"
Ground to Top of Seat	29" - 31"
Ground to Top of Deck	8.75"
Deck to Top of Seat	19.75" - 21.75"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	40"
Weight of Unit: (assembled)	302 lbs
Base	152 lbs
Seat w/Arms	44 lbs
Batteries	52 lbs ea
Tires:	Pneumatic
Front	13"
Rear	13"
Standard Seat:	Captain's Pillow Top High-Back
Slide Seat	Yes
Back Height - no headrest	20"
Back Height with headrest	26.5" - 30"
Width x Depth	20" x 18"
Color	Black (coming soon)
Warranty:	
Frame	5 Years
Drive Train	1 Year
Electronics ²	1 Year
Batteries	1 Year

Battery range at 200 lbs., and will vary due to rider weight, drive surface and terrain. Batteries not included.



New button covers, hand (parking) brake lock, LED battery meter.



Updated control panel on wraparound



New 13" pneumatic tires.



New ultrabright LED rear lights.



² Electronics warranty excludes batteries.