

# HERITAGE COLLECTION

## COMFORT AND STYLE

FDA Class II  
Medical Device<sup>1</sup>

### Superior Design, Exceptional Strength

Comfort and style complement one another perfectly in the Heritage Collection from Pride® Power Lift Recliners. With a superior design that provides exceptional strength, Heritage Collection delivers style and durability that will last for years to come. Pride's Heritage Collection offers numerous standard features and quiet and smooth performance.



Footrest Extension

USB Phone Charger



Scan here for more information on the **Heritage Collection.**



FOLLOW US ON SOCIAL MEDIA

**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®

### FEATURES & BENEFITS

- Easily adjustable leg levelers provides additional support for uneven surfaces
- Easy-to-use hand control with USB charger port (not available when adding optional heat and massage)
- Large dual pockets for extra storage
- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (except LC-358XL and LC-358XXL)
- Quiet and smooth lift and scissor mechanism system
- Integrated battery backup ensures chair is operational when power goes out
- Energy-efficient power box with 9' cord
- Premium furniture grade laminate/hardwood frame for added strength and durability

### OPTIONS\*

- Heat and massage (except XXL)
- Cleaning kit
- Head, lumbar and bolster pillows
- First Up<sup>®</sup> technology (Not available on LC-358XL & LC-358XXL)
- Non-skid leg levelers
- Lithium battery backup

\* Additional charges apply for all options. Extra lead time may apply for some options.

MODEL	LC-358PW <sup>1,4</sup>	LC-358S <sup>1,4</sup>	LC-358M <sup>1,4</sup>	LC-358L <sup>1,4</sup>	LC-358XL <sup>1,3,5</sup>	LC-358XXL <sup>1,3,5,6</sup>
POSITION	3-Position 5'3" & Below	3-Position 5'3" & Below	3-Position 5'4" - 5'9"	3-Position 5'9" - 6'0"	3-Position 5'8" - 6'1"	2-Position 5'8" - 6'1"
SIZE	Petite Wide	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
SEAT-TO-FLOOR	19"	19"	20"	21"	20"	20"
SEAT DEPTH	19"	19"	20"	22"	22"	22"
SEAT WIDTH	24"	18"	20"	22"	26"	30"
TOP OF BACK TO SEAT	23.5"	24"	24.5"	25.5"	29"	29"
BACK STYLE	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes	Yes
HEAT AND MASSAGE	Optional	Optional	Optional	Optional	Optional	N/A
SEATING TYPE	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating
FOOTREST EXTENSION	Yes	N/A	Yes	Yes	N/A	N/A
DISTANCE FROM WALL	18"	18"	17"	19"	19"	14"
OVERALL WIDTH	40"	35"	35"	37"	39.5"	43"
CHAIR WEIGHT	130 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
WARRANTY	7-Year Pro-rated <sup>7</sup>	7-Year Pro-rated <sup>7</sup>	7-Year Pro-rated <sup>7</sup>	7-Year Pro-rated <sup>7</sup>	7-Year Pro-rated <sup>7</sup>	7-Year Pro-rated



Cloud 9 - Quick Ship

Lexis Sta-Kleen<sup>®</sup>  
OPTIONAL

UltraFabrics<sup>®</sup> Ultraleather<sup>®</sup>  
OPTIONAL

Crypton<sup>®</sup> Aria<sup>2</sup>  
OPTIONAL

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

1. Pride<sup>®</sup> Power Lift Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.  
 2. Crypton<sup>®</sup> Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.  
 3. Removable, fully padded chaise.  
 4. Unique spring design and a one-piece soft foam seat for outstanding relaxation.  
 5. Features Dual Motor Sync System.  
 6. Only available in Cloud 9 Stone, Cloud 9 Walnut, Lexis Sta-Kleen<sup>®</sup> Chestnut and Lexis Sta-Kleen<sup>®</sup> Mushroom.  
 7. Lifetime warranty on steel lift frame.

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.