

The Essential Collection from Pride® Power Lift Recliners provides exceptional value and outstanding comfort. Several fabric options will match any décor while delivering rest and worry-free operation.



Easy-to-use hand control

ENERGY EFFICIENT WITH BUILT-IN BATTERY BACKUP



1) Pride® Power Lift Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Scan here for more information on the **Essential Collection.**



FOLLOW US ON SOCIAL MEDIA

## FEATURES & BENEFITS

- Easy-to-use hand control adjusts your chair to find your perfect position
- Large dual pockets offer ample storage space that is easily accessible
- Smooth and durable lift and scissor mechanism for long-lasting and quiet performance
- Removable back lets you adjust chair padding as needed
- Integrated battery backup ensures chair is operational when the power goes out
- Energy-efficient power box with 9' cord
- Lifetime warranty on steel lift frame
- Fire retardant (California specifications; bulletin 117-2013)

## OPTIONS\*

- Head, bolster, and lumbar pillows (LC-250 only)
- Cleaning kit
- Non-skid leg levelers

\* Additional charges apply for all options.

MODEL	LC-250 <sup>1</sup>	LC-105 <sup>1</sup>
POSITION	3-Position 5'4" - 5'9"	3-Position 5'3" and below
SIZE	Medium	Medium
WEIGHT CAPACITY	375 lbs.	325 lbs.
SEAT-TO-FLOOR	20"	18"
SEAT DEPTH	20"	19"
SEAT WIDTH	22"	20"
TOP OF BACK TO SEAT	26"	28"
FULLY PADDED CHAISE	Yes	Yes
SEATING TYPE	Total Comfort	Sinuous Spring
DISTANCE FROM WALL	20"	17"
OVERALL WIDTH	34"	33"
CHAIR WEIGHT	118 lbs.	96 lbs.
WARRANTY	2-Year Limited	1-Year Limited



LC-250<sup>1</sup>



LC-250 FABRICS **QUICK SHIP**

LC-105 FABRICS **QUICK SHIP**

(US) 800-800-8586 / (CAN) 888-570-1113 / [pridemobility.com](http://pridemobility.com)

1) Pride® Power Lift Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>