Jazzy® Carbon¹

K0899 Power Mobility Device, not coded by DME PDAC or does not meet criteria
300 lbs.
24"
24"
37"
Up to 3.7 mph
2"
7"
8.5" solid
N/A
24v/250W brushless motors
Regenerative and electromechanical
Front
17" x 17" black foam
20.25" with standard cushion 18.75" without cushion
Up to 9.3 miles
35A Brushless Micon
Base: 39.6 lbs. Battery: 4 lbs.
43.6 lbs. with batteries
24V 12AH (288WH) lithium-ion
2 amp, lithium-ion off-board
Lifetime limited on frame; 1-year limited on electrical control/ lever system, motor/drive system, bearing sleeve and battery

- 1 Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).
- 4 Total weight includes base, seat and batteries.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.






(US) 800.800.8586 • (Canada) 888.570.1113 pridemobility.com | • • •

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or -3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



AMERICA'S FAVORITE POWER CHAIR pridemobility.com



UNFOLD, SNAP, AND GO!

This is our most portable Jazzy® yet. In just seconds, the Jazzy Carbon folds flat for easy transport, and unfolds again to get you moving. Ultra-light but exceptionally strong, this electric wheelchair is your ultimate travel companion.

Weighing in at a stunning 43.6 pounds, the Jazzy Carbon still provides top speeds up to 3.7 mph with a weight capacity of 300 pounds.

Features

- Removable joystick with built in USB charger
- 24V 12Ah (288WH) airline compliant lithium battery*
- Front suspension for a smooth ride
- Ultra-lightweight, carbon fiber frame specifically designed for traveling
- Compact base easily folds and unfolds for transport in a few simple steps
- Padded seat back with storage pocket and soft, removable seat cushion for all-day comfort
- Free-wheel levers easily disengage the brake, allowing this lightweight power chair to be pushed like a manual chair
- Under-seat storage and adjustable lap belt

Available Colors





Black

White

Accessories

- Stealth® Simplicity cushion
- Military patches
- Additional battery

* Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® Carbon and its lithium ion battery are permitted.

