

RevVA3

REGAL

FULL-SIZE



merits

3/4/2020

FEATURES

- 10" knobby drive tires for superior traction
- 3" ground clearance for easy obstacle clearance
- Dynamic LiNX controller for easy programming and diagnostics
- 16", 18" or 20" ergonomic, semi-reclining seats: up to 135°

SPECIFICATIONS

Weight Capacity	300 lbs
Max Speed	5 mph
Per-Charge Range	18 miles
Batteries	U1-34Ah x 2
Ground Clearance	3"
Turning Radius	21"
Overall Length	39"
Overall Width	24"
Drive Wheels	10" foam filled tires
Rear Casters	8" foam filled tires
Seat-to-Floor Height	20"-22" (captain), 16.5"-18.5" (pan)
Seat Width	16", 18" or 20"
Seat Depth	16" or 18"
Back Height	18" or 20"
Controller	Dynamic LiNX
Warranty	18 Months : Motor & Electronics
	5 Year/ Limited Lifetime Frame
	1 Year/ Batteries

OPTIONS

- Swing-away joystick
- USB charger
- Oxygen tank holder
- Cane/Crutch holder
- Walker holder
- Foldable rear basket
- Weather cover
- Elevating legrests

MODELS

Regal w Captain Seat	P310-2S C
Regal w Pan Seat	P310-2S S



Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-32

merits