

Roadster



S731RS3-SBB / Galaxy Blue
S731RS3-ARB / Cherry Red

SPECIFICATIONS

| | |
|---------------------------|----------------------------|
| Product # | S731RS3 |
| Weight Capacity | 300 lbs. |
| Range | up to 12 mi. |
| Seat Width | 18" |
| Seat Depth | 18" |
| Backrest Height | 14.5" from seat |
| Overall Length | 40.5" |
| Overall Width | 25.5" |
| Overall Height | 36.6" |
| Turning Radius | 47" |
| Ground Clearance | 2.2" |
| Total Weight | 103 lbs. |
| Weight of Heaviest Piece | 32 lbs. |
| Max Speed | 4 mph |
| Basket Size | 12.5wx10hx6d |
| Motor | 24V 180W |
| Battery Size (2 required) | 12V/12AH |
| Charger | 29V2AH |
| Gradient | 6° |
| Front Wheel (s) | 7" |
| Rear Wheels | 8" |
| HCPCS Code | na |
| Standard Colors | Cherry Red and Galaxy Blue |



- Maximum user weight 300lbs.
- Maximum speed 4 mph.
- Drive range up to 12 miles.
- 3 wheel design for easy maneuverability.
- Disassembles into five pieces for easy transportation.
- Weight of heaviest piece - 32 lbs.
- Headlight to help visibility.
- Convenient storage pocket.

MSRP **\$2,090**



4245 Evans Ave., Fort Myers, FL 33901

239-772-0579

W www.meritsusa.com E sales@meritsusa.com