



6Runner 10

The 6Runner 10 is a mid-size power chair that is designed for the daily user. It sports a mid-wheel drive system and has the added stability of having six points of contact with the ground. The 6Runner 10 also has our stable four post seating system and suspension system for smooth drive operation. The unit has easy battery access, an 18in. wide reclining seat and 300 lb. weight capacity.

- + 300lb Weight Capacity
- + Mid-wheel Drive
- + Full Suspension
- + Easy Battery Access
- + Four Post Seating System



Model Number: 888WNLM
Medicare Billing Code: K0823

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

SPECIFICATIONS

WEIGHT CAPACITY	300 lbs.
WIDTH ¹	23"
LENGTH ¹	42"
UNIT WEIGHT ³	191 lbs.
MAX SPEED	5 mph

GROUND CLEARANCE	3"
SAFE CLIMBING ANGLE	6° 10%
FRONT REAR WHEELS	7" Solid 7" Solid
DRIVE WHEEL	10" Solid
SUSPENSION TYPE	Limited

BATTERY OPTIONS	U1
BATTERY WEIGHTS	21.3 lbs. x 2
RANGE ² UP TO	10 miles
CHARGER	5A Off-board
STFH	21.5" - 23.5"

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.