



6Runner 14

The 6Runner 14 is a heavy duty power chair that is designed for the daily user. It sports a mid-wheel drive system and has the added stability of having six points of contact with the ground. The 6Runner 14 also has our stable four post seating system and suspension system for smooth drive operation that will handle most driving surfaces. The unit has a 450 lb. weight capacity, reclining seat, adjustable arm rest and easy battery access.

- + 450lb Weight Capacity
- + Mid-wheel Drive
- + Four Post Seating System
- + Quick Battery Access
- + Reclining Captain Seat

● Burgundy

Model Number: 888WNLLHD
Medicare Billing Code: K0825

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	42" x 26" x 49"
Number, Size of Tires	Front, Rear	2, 7" (2, 14") 2, 7"
Suggested User Weight ¹	lb	Max. 450
Weight of Heaviest Part ²	lb	-
Weight ³	lb	277
Each Battery Module Weight	lb	34
Battery Size	Type x PCS	12V50AH x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	1.34 x 2
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	10(6)
Range (Per Charge) ⁶	Mile	21.3 - 23.8
Turning Radius	in	18"
Model #		888WNLLHD
Standard Seat	W x D x BH	SV2121211/2RB*

**Additional seat sizes and options available as an upgrade. Please inquire for more details.*

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.