

# Jimmie

The Jimmie is a portable and light weight power chair designed for the user on the go. The unit has a quick-disconnect system for effortless disassembly, advanced electronics and high efficiency in-line motors. The Jimmie disassembles effortlessly into three pieces, with the heaviest weighing just 35.8 lb. It also has all the features of larger power chairs including height & angle adjustable arm rests, adjustable footplate and adjustable seat height.



- + 250lb Weight Capacity
- + Wireless Take-A-Part technology
- + Adjustable Seat Height
- + Easy to Dismantle into Three Parts
- + Padded Arm Rests



Burgundy

Model Number: UL8WPBS

Medicare Billing Code: K0821

HCPSC codes provided should not be considered as legal advice, and does not guarantee reimbursement.

## SPECIFICATIONS

WEIGHT CAPACITY	250 lbs.
WIDTH <sup>1</sup>	22.5"
LENGTH <sup>1</sup>	35"
UNIT WEIGHT <sup>3</sup>	106 lbs.
MAX SPEED	3.75 mph

GROUND CLEARANCE	2"
SAFE CLIMBING ANGLE	6°   10%
FRONT   REAR WHEELS	6" Solid
DRIVE WHEEL	9" Solid
SUSPENSION TYPE	Limited

BATTERY OPTIONS	12V x2 Batt. Pack
BATTERY WEIGHTS	21.3 lbs.
RANGE <sup>2</sup> UP TO	10 miles
CHARGER	2A Off-board
STFH	19.5" - 20"

**Note:** <sup>1</sup> On level Road. <sup>2</sup> Without basket and seating. <sup>3</sup> With basket, seat & batteries. <sup>4</sup> Base on user weight of 200LB. <sup>5</sup> Based on user weight of 165LB. <sup>6</sup> After the batteries and mechanical moving parts fully break in.