Remote Operated

Transformer Scooter

- No Tools Required
- Folds Compact & Portable
- Up to 13.5 Miles Per Charge
- Lightweight
- Easy to Carry
- 300lb Weight Capacity
- Remote Operated
- Ideal for Travel
- Fits in Most Vehicles
- Airline Approved Battery

Shown in Yellow Available in Four Colors

RANSFORMER FOLDING SCOOTER





| Specifications | Transformer Scooter |
|---------------------|--------------------------|
| Overall Length | 37.4″ |
| Overall Width | 22.25" w/ Armrests |
| Tires | Front – 6x1.5" |
| | Rear – 7x2.36" |
| Maximum Speed | 3.7mph |
| Range | Up to 13.5 miles |
| Turning Circle | 55″ |
| Weight (w/ battery) | 53lbs |
| Color | Red, Blue, Yellow, Black |
| Brake System | Intelligent Regenerative |
| | Electromagnetic Braking |
| Drive System | Rear Wheel Drive |
| Handlebar | PU Foam |
| Weight Capacity | 300lbs |
| Seat Width | 17" |
| Seat Height (from | 21.75″ |
| ground) | |
| Battery | 24V10Ah Lithium |
| Motor | 120W |
| Battery Charger | DC24V/2A |
| Folded Dimensions | 17" x 22.3" x 25" |