



## Streamer Sport

The Streamer Sport is a mid-size power chair designed for the daily user. It has a 300 lb. weight capacity and rear-wheel drive maneuverability. The unit comes standard with height and width adjustable arm rests, solid tires, adjustable angle footplate and a 19in. captains seat. The Streamer Sport has a 10-mile battery range and 3in. ground clearance that will tackle daily use with ease.

- + 300lb Weight Capacity
- + 10" Drive Wheels
- + Adjustable Angle Foot Plate
- + Height and Width Adjustable Arm Rests
- + Captain Seat

 Red

Model Number: 888WA  
Medicare Billing Code: K0823

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

## Technical Specifications

Overall Dimensions (L x W x H)	in	36" x 24" x 45"
Number, Size of Tires	Front, Rear	2, 6" (2, 10")
Suggested User Weight <sup>1</sup>	lb	Max. 300
Weight of Heaviest Part <sup>2</sup>	lb	-
Weight <sup>3</sup>	lb	203
Each Battery Module Weight	lb	23
Battery Size	Type x PCS	U1 x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	.60 x 2
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle <sup>4</sup>	%(°)	10(6)
Range (Per Charge) <sup>6</sup>	Mile	10
Turning Radius	in	17.7"
Model #		888WA
Standard Seat	W x D x BH	SV191821B

**Note:** <sup>1</sup> On level Road. <sup>2</sup> Without basket and seating. <sup>3</sup> With basket, seat & batteries. <sup>4</sup> Base on user weight of 200LB. <sup>5</sup> Based on user weight of 165LB. <sup>6</sup> After the batteries and mechanical moving parts fully break in.