

TRANSFORMER 2


TRANSFORMER 2 FOLDING SCOOTER

Remote
Operated

- Contemporary New Look
- Folds Compact & Portable
- Up to 13.5 Miles Per Charge
- Lightweight
- Easy to Carry
- 300lb Weight Capacity
- Remote Operated
- Ideal for Travel
- Fits in Most Vehicles
- Airline Approved Battery
- Front & Rear Suspension
- LED Headlight



Shown in Blue
Available in Four Colors



| Specifications | Transformer S3026 |
|---------------------------|---|
| Overall Length | 38" |
| Overall Width | 18" w/ Armrests |
| Tires | Front – 6x1.5" Rear – 7x2.36" |
| Maximum Speed | 3.75mph |
| Range | Up to 13.5 miles |
| Turning Circle | 55" |
| Weight (w/ battery) | 55lbs |
| Color | Red, Blue, Yellow, Black |
| Brake System | Intelligent Regenerative Electromagnetic Braking |
| Drive System | Rear Wheel Drive |
| Ground Clearance | 3 1/4 " |
| Weight Capacity | Up to 300lbs |
| Seat Width | 17" |
| Seat Height (from ground) | 22" |
| Battery | 24V10Ah Lithium Ion |
| Motor | 120W |
| Battery Charger | DC24V/2A |
| Folded Dimensions | 18" x 22" x 25" |