



Enhance Your Life

TRIAXE Cruze

Why settle? When you can have it all!

- Convenience & Portability of a Folding, Travel Scooter
 - Comfort & Performance of a Full-Size Scooter
 - Folds Flat for Transport & Storage
- Durable Aluminum Construction – 350lb Weight Capacity
- Electromagnetic Brakes, Free Wheel Mode & Speed Control
 - Deluxe Basket Included
- Speeds to 6mph & Up to 18 Miles per Charge



Model	TRIAXE Cruze T3055
Overall Length	39"
Overall Width	19.5"
Tires	8" Solid Front Tire 9.5 Solid Rear Tires
Maximum Speed	Up to 6 mph
Range	Up to 18 miles
Turning Radius	30"
Weight (without battery pack)	58lbs
Color	Fiery Red, Reflex Blue
Brake	Electric Brake Systems
Drive System	Rear Wheel Drive, Brushed Motor
Handlebar	Aluminum with Rubber Hand Grips
Maximum Weight Capacity	Up to 350lbs
Seat Width	18" standard
Seat Height	23"
Battery	24V 12AH
Motor	24v 180w
Battery Charger	24V 2AH