

Bobcat X3

HCPCS:




FIGURE A




FIGURE B




FIGURE C



FIGURE D




FIGURE E



FIGURE F

- Ultra-lightweight for superior portability
- Quick and easy disassembly
- Color-through shroud hides scuffs and scratches
- Angle adjustable flip-up armrests for maximum comfort and convenient boarding and exiting
- Easy finger-pull throttle control
- Height adjustable seat with fold down backrest
- Flat-free, non-marking tires for worry-free travel
- Large, metal carry basket
- Easy-to-adjust tiller angel places the scooter's controls in the ideal position
- Anti-tip wheels for additional safety
- Easy free wheel operation

PRODUCTS

Description

BOBCATX3	BOBCAT X3 3-WHEEL SCOOTER, RED, 1/EA
----------	--------------------------------------

SPECIFICATIONS

Description

Base Weight	50 lbs
Batteries (pair)	12V x 10AH
Battery Weight (pair)	18 lbs
Brakes	Electromagnetic
Charger	1.5A Offboard
Climbing Angle	6°
Controller	45A
Dimensions	42" x 19.25"
Floor to Seat Height	17" - 19"
Ground Clearance	2.5"
Heaviest Piece Weight	28 lbs
Maximum Range	7.5 miles
Motor	24V x 270W x 4700 rpm
Seat Dimensions	16.5" x 13.5"
Seat Weight	16 lbs
Top Speed	4 mph
Total Weight (excl. hand control)	84 lbs
Turning Radius	45"
Warranty on Batteries	12 months
Warranty on Drive Train Components	24 Months

SPECIFICATIONS	Description
Warranty on Electronic Controller	24 Months
Warranty on Frame	Lifetime
Weight Capacity	265 lbs
Wheels	8" x 2" Flat Free
ACCESSORIES	Description
SF8050	Sun Shade for Scooter,1/ea, CS
SF8020	Rear Basket Set, LType, 1/ea, CS
SF8010	Oxygen Holder Set, L Type,1/ea, CS
SF8030	Crutch Holder Set, SType, 1/ea, CS
STDS1040S	Cup Holder, Swivel,1BX,6/CS, CS
PARTS	Description
S31003	Grip,Phoenix,Odyssey,Day-3/4,1/ea, CS
S350016910	Rear Axle Bolt,Phoenix, Bobcat X, 1/ea, CS
S31026	Anti Tip Wheel,Phoenix/Bobcat,1/ea, CS
S31052	Armrest Knob POV 1/ea, CS
BC35051	BasketHook Set,Phoenix,Scout,Ventura,1/e, CS
BC40AFUSE	Bobcat Fuse 40A,1ea, CS
BC31132	Rear Tire, Black Rim, Bobcat , 1/ea, CS
BC31142	Front Wheel ,Black Rim, Bobcat 3, 1/ea, CS
BC31141	FrontAxleAssemblyBobcat3,1ea, CS
C12-071-00100	Speed Dial, Bobcat X3 / X4, 1/ea, CS
C08-050-00301	Upper Tiller Cover, Bobcat X3 / X4, 1/ea, CS
C18-071-00100	Key Switch, Bobcat X3 / X4, 1/ea, CS
C18-071-00200	Key, Bobcat X3 / X4, 1/ea, CS
C08-071-00100	Lower Tiller Cover, Bobcat X3 / X4, 1/ea, CS
C08-069-01001	
C02-071-00100	Tiller, Bobcat X3/ X4, 1/ea, CS
C02-069-00401	TillerSwivelAssembly,Bobcat X3/ X4, 1/ea, CS
C29-016-00301	TillerBearingSet,BobcatX3/X4, 1 set/pk, CS
C08-071-00400	Rear Shroud, Bobcat X3 / X4, 1/ea, CS
C08-071-00200	Main Shroud, Bobcat X3 / X4, 1/ea, CS