

SHOPRIDER®



K0806

DASHER 9 [GK-73]

The Dasher® 9 is fun for traveling! The Dasher® 9 is equipped with "Splash-n-Dash" quick color change moldings to give you the option of changing the color of your scooter effortlessly. The versatile Dasher® 9 is stylish, reliable, and provides a comfortable ride wherever you want to go!

Toll-Free: 800.743.0772
www.shoprider.com



Adjustable Tiller

MID SIZE SCOOTER

- 300 lb. weight capacity
- Quick connect battery pack
- Charge battery pack at any convenient location on/off the scooter
- Puncture proof tires
- Adjustable seat height
- Front Basket included

DASHER 9



Overall Dimensions (L x W x H)	in	38.5 x 22.5 x 36
Number, Size of Tires	Front (rear)	1, 9" (2, 9")
Suggested User Weight (on the Level Road)	lb	Max. 300
Weight of Heaviest Part (Front Chassis Assembly)	lb	35
Total Weight (w/ 21Ah Battery Pack, Basket, Seat)	lb	110
Battery Pack Module Weight	lb	35
Battery Module Capacity	_V_Ah x Pcs	12V21Ah x 2
Charger		Off-Board (2A)
Output Power of Motor (Reference Only)	hp	0.75
Maximum Speed ²	mph	3.8
Range ⁴ (per charge w/ STD 21AH battery /200 lbs User Weight) (After the battery and mechanical moving parts fully break in)	mile	11.2
Turning Radius	in	31.5
Handle Bar	Type	T-bar

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.