

Explorer

Lightweight Folding Scooter

eFOLDi 

Distributed in the USA by  **AFIKIM**
ELECTRIC VEHICLES

Features & Benefits

- Lightweight folding scooter using a moulded aluminium frame.
- Folds and unfolds in seconds and needs no assembly.
- Weighs just 38lbs (without battery) making it easy to transport.
- Compact size means it can be stowed almost anywhere and fits the trunk of any car.
- Low centre of gravity and pneumatic tyres give a safe, smooth and stable ride.
- Airline approved battery for ease of travel.
- Tight turning circle for easy indoor use (eg grocery store).
- Powerful 240W motor allows use on more challenging terrain.



Specifications

Size & Weight

Folded Size	26 x 22 x 15.6 inches
Unfolded Size	47.6 x 22 x 32.6 inches
Seat Height	19 inches
Footrest Height	8.5 inches
Ground Clearance	4.7 inches
Front Wheel Size	12 inches
Rear Wheel Size	10 inches
Weight (without battery)	38lbs
Weight (with battery)	42lbs

Technical Specifications

Motor	240W 24VDC Hub Motor
Battery	24V 12Ah Li-Ion, 288 Wh
Charger Input	Input 110-220VAC, 50-60Hz
Charger Output	29.4V 2A

Performance

Use	Pavement, Indoors & Roads
Speed	4mph / 8mph
Maximum Carry Capacity	265lbs
Maximum Incline	1:6 Gradient (9.5 degrees)
Turning Radius (Left/Right)	40 inches (L) – 33 inches (R)
Maximum Kerb Height	1.1 inches
Maximum Range	14 Miles

Brake System	Mechanical handbrake
Throttle	Twist grip
Frame Material	Aluminium Alloy
Lights	Headlamp, taillights, brake & turn indicators

For a full range of support resources please visit https://efoldi.com/us_explorer/