

Ultra Lightweight Folding Mobility Scooter

Distributed in the USA by  **AFIKIM**
ELECTRIC VEHICLES

Features & Benefits

- One of the world's lightest folding mobility scooters using a revolutionary magnesium alloy frame.
- Folds and unfolds in seconds and needs no assembly.
- Weighs just 33lbs (without battery) making it easy to transport.
- Compact size means it can be stowed almost anywhere and fits the trunk of any car.
- Low centre of gravity and pneumatic tyres give a safe, smooth and stable ride.
- Airline approved battery for ease of travel.
- Tight turning circle for easy indoor use (eg grocery store).
- Fully FDA approved.



Specifications

Size & Weight

Folded Size	27.5 x 22 x 15 inches
Unfolded Size	46 x 22 x 37 inches
Seat Height	20 inches
Footrest Height	4.7 inches
Ground Clearance	2.7 inches
Front Wheel Size	8 inches
Rear Wheel Size	10 inches
Weight (without battery)	33lbs
Weight (with battery)	38lbs

Performance

Use	Pavement and Indoors
Speed	4mph
Maximum Carry Capacity	265lbs
Maximum Incline	1:10 Gradient (6 degrees)
Turning Radius	40 inches
Maximum Kerb Height	1.5 inches
Maximum Range	10 Miles

Technical Specifications

Motor	180W 24VDC Hub Motor	Brake System	Electromagnetic
Battery	24V 10aH Li-Ion, 240 Wh	Throttle	Thumb paddle
Charger Input	Input 110-220VAC, 50-60Hz	Frame Material	Magnesium Alloy
Charger Output	29.4V 2A	Lights	None