



4 wheel GA541

Seat Color Charcoal Vinyl

Body Colors Silver Re

The Golden Avenger four-wheel, heavy duty scooter provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of nature's most challenging terrains. Featuring a 500 lb. weight capacity.





GA-541

Avenger 4 wheel **Specifications**

Model	GA541
Medicare Code	Not Classified
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	7 mph
Operating Range 1	18 miles
Ground Clearance to Center D	
Turning Radius	67"
Type Batteries	2-NF22
Colors	Silver, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PGS 140 Amp
Dimensions:	
Length	57.5"
Width	24.5"
Height (ground to top	
of back of seat)	42.5" - 45.5"
Height (ground to top	
of back of headrest)	50" - 56"
Ground to Top of Seat	25" - 28"
Ground to Top of Deck	6.5"
Deck to Top of Seat	18.75" - 21.75"
Number of Seat Height Posit	ions 4
Size of Increments	1"
Front Axle to Rear Axle	38.5"
Weight of Unit: (assembled) 268 lbs
Front Half	68 lbs
Rear Half	68 lbs
Base	N/A
Seat w/Arms	55 lbs
Batteries	37 lbs ea
Tires:	Pneumatic
Front	13"
Rear	13"
Casters	N/A
Rear Anti-Tippers	3"
Standard Seat:	Captain's Seat
Slide Seat	N/A
Back Height - no headrest	21" - 24"
Back Height with headrest	27.5" - 31.5"
Width x Depth	20" X 17"
Color	Gray Vinyl
Warranty:	dray viriyi
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Batteries	1 Year

Battery range will vary due to rider weight, drive surface and terrain. Foam-filled tires are optional. The operational mileage will be less if scooter is equipped with foam filled tires. Batteries not included.
² Electronics warranty excludes batteries.



Features

- 500 lb weight capacity
- Speeds up to 7 mph
- Bright safety lights
- 22 NF batteries
- Up to 18 Mile Range
- 4-Way Flashers Great Warranty
- Oversized tires
- Removable basket





