



Smart Fold, Easy Go

Luggie Standard

LENGTH	39 in (99.1 cm)
WIDTH	18 in (45.7 cm)
LENGTH (FOLDED)	26 in (66 cm)
WEIGHT CAPACITY	250 lbs. (113.4 kg)
WEIGHT (SCOOTER)	51.7 lbs (23.5 kg)
WEIGHT (SCOOTER+BATTERY)	56.9 lbs (25.8 kg)
BATTERY	1x 8.5AH Lithium Ion
BATTERY CHARGER	5-amp off-board charger
MAX. SPEED	4 mph (6.4 kph)
RANGE*	11 miles (17.7 km)
TURNING RADIUS	35.5 in (90.2 cm)
GROUND CLEARANCE	2.5 in (6.4 cm)
MAX. INCLINE	6 degrees
WHEELS/TIRES	Front: 2" x 6" PU solid Rear: 2" x 7" PU solid
DRIVE SYSTEM	24-volt DC motor, sealed mini-transaxle, rear-wheel drive

STANDARD COLORS



The information contained herein is correct at the time of publication. Factors such as battery charge, battery condition, user weight, tire pressure, and terrain type will cause variation to maximum range and speed. Due to continual product improvement and manufacturing tolerances; length, width, turning radius and ground clearance can be subject to variance of +/- 3%. We reserve the right to alter specifications without prior notice.



Lightweight • Easy-folding • Travel-friendly

The Original Travel Mobility Scooter. Innovative folding design does not require tools. Lithium ion battery is light-weight and airline approved.

