



Item #	Seat	UOM
PROWLER3310MG20CS	20" Captain's Seat	1/cs
PROWLER3310MG22CS	22" Captain's Seat	1/cs



HCPCS
K0899

Prowler 3310

Large Scooter

- Armrests are padded and adjustable.
- Ergonomic throttle control.
- Height-adjustable swivel seat features fold-down backrest.
- Seat back reclines and features adjustable headrest.
- Available with 20" or 22" width Captain's Seat
- Seat position can be adjusted forward or backward.
- Pneumatic tires and full suspension.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Programmable controller.
- Anti-tip wheels.
- Headlight.
- Rearview mirrors.
- Turn signals/hazard lights.
- Delta tiller.
- Easy freewheel operation.
- Additional hand-operated drum brake.

specifications

HCPCS	K0899	Batteries (pair)	12V x 50AH
Weight Capacity	500 lbs.	Charger	5A Offboard
Top Speed	7 mph	Brakes	Electromagnetic & Drum
Maximum Range	23 miles	Front Wheels	12" x 3.5" Pneumatic
Turning Radius	50"	Rear Wheels	12" x 3.5" Pneumatic
Climbing Angle	7.5°	Base Weight	130 lbs.
Ground Clearance	5"	Battery Weight	66 lbs.
Dimensions	26" (W) x 51" (L)	Seat Weight	55 lbs.
Floor-to-Seat Height	25" - 28"	Total Weight	251 lbs.
Seat Type	Captain's	Warranty on Frame	Lifetime
Seat Dimensions	20" (W) x 20" (D) 22" (W) x 20" (D)	Warranty on Electronic Controller/Drive Train Components	14 months
Motor	24V x 650W x 4350 rpm	Warranty on Battery	6 months