

# Super Micro Electric Scooter

## S539 Mini-Coupe

### Technical Specs

max speed : up to 4.5mph  
 range : up to 10 mi.  
 turning radius : 34"  
 ground clearance : 3"  
 motor : DC 24V, 110W  
 controller : S-Drive 70A  
 battery : 12V / 17AH x 2pcs  
 charger : 1.5A off-board  
 gradient : 6°  
 front wheel : 8" foam filled  
 rear wheel : 9" foam filled  
 brake : intelligent, regenerative, electromagnetic brakes

### Features

- Auto latching lockup mechanism for simple frame separation using only one hand.
- All new frame design easily disassembles for convenient transport and storage.
- Convenient 2 x 18Amp battery pack provides increased driving range.
- 300-lb weight capacity.
- Color available : Red.

### Weight (lb)

base	battery/ea	seat	tl wt w/ batteries	capacity
53	11	15.4	90	300

### Dimensions (in)

length	width	height	wheel base	seat width
39.4	20	32	30	16.5
seat depth	back height	seat height from deck	seat height from ground	
15.5	14	16.5	18.5	



Medicare Code  
K0800

