

# Echo Folding

The Echo Folding is designed for the user on the go. The scooter folds for transportation and unfolds for use in a matter of seconds. The unit weighs in at just 59 lb. and can store easily in small spaces. It comes standard with solid tires, electric brake system and a curved tiller for added comfort. Available in red, gray, black and yellow, the Echo Folding will get you to your destination in style.



- + 200lb Weight Capacity
- + Adjustable Tiller
- + Solid Tires
- + 59 lb Overall Weight



- Red
- Black
- Yellow

**Model Number: FS777**

HCPSC codes provided should not be considered as legal advice, and does not guarantee reimbursement.

## SPECIFICATIONS

WEIGHT CAPACITY	200 lbs.
WIDTH	17.5"
LENGTH	39.5"
UNIT WEIGHT <sup>1</sup>	59 lbs.
MAX SPEED <sup>2</sup>	3.75 mph

GROUND CLEARANCE	2"
SAFE CLIMBING ANGLE	6°   10%
FRONT   REAR WHEELS	6" Solid
DRIVE WHEEL	8" Solid
SUSPENSION TYPE	Limited

BATTERY OPTIONS	12V10Ah
BATTERY WEIGHTS	7 lbs.
RANGE UP TO <sup>3</sup>	9 miles
CHARGER	2A Off-board
TILLER TYPE	T-Bar

Note: <sup>1</sup> Maximum unit weight dependent on seating configuration and battery choice, standard equipment weight shown. <sup>2</sup> Maximum speed may vary, user weight and road grade may affect top speed. <sup>3</sup> Range can vary depending on level of charge, user weight and road conditions. HCPSC codes provided should not be considered as legal advice, and does not guarantee reimbursement.