

Ventura 4

HCPCS: K0807-1, K0807



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Seat position can be adjusted forward or backward
- Front and rear suspension
- Programmable controller
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Delta tiller
- Easy free wheel operation

PRODUCTS

PRODUCTS	Description
VENTURA418FS	VENTURA 4-WHEEL 18" FOLDING SEAT, RED/BL

SPECIFICATIONS

SPECIFICATIONS	Description
Base Weight	99 lbs.
Batteries (pair)	12V x 33AH
Battery Weight (pair)	28 lbs.
Brakes	Electromagnetic
Charger	4A Offboard
Climbing Angle	7.5°
Dimensions	48" (L) x 25" (W)
Floor to Seat Height	22" - 24"
Front Wheels	10" x 3" Flat-free
Ground Clearance	2.5"
HCPCS	K0807
Heaviest Piece Weight	70 lbs.
Maximum Range	16 miles
Motor	24V x 350W x 5100 rpm
Rear Wheels	10" x 3" Flat-free

SPECIFICATIONS	Description
Seat Dimensions	18" (W) x 17.5" (D)
Seat Type	Folding
Seat Weight	26 lbs.
Top Speed	6 mph
Total Weight (excl. hand control)	153 lbs.
Turning Radius	54"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Frame	Lifetime
Weight Capacity	400 lbs
Weight Capacity (13046 13046-HD)	400 lbs.

ACCESSORIES	Description
SF8030	Crutch Holder Set, SType, 1/ea, CS
SF8010	Oxygen Holder Set, L Type,1/ea, CS
SF8050	Sun Shade for Scooter,1/ea, CS
SF8020	Rear Basket Set, LType, 1/ea, CS

PARTS	Description
VENT-86	Shroud Knob, Ventura, 1/ea, CS
LRM412123	21AH Battery,POV, 1/ea, CS
VENT-15	ShockAbsorberVentura3&4, CS
VENT-14	FrontBumperAssemblyVentura4,1/ea, CS
VENT-44	Head Light Assembly, Ventura 4, 1/ea, CS
VENT-17	FrontFrameVentura4,1/ea, CS
VENT-16	Anti-Tipper Wheel, Ventura 3 & 4,1/ea, CS
VENT-11	MiddleShroudVentura3&4,1/ea, CS
VENT-09	BatteryStrapsVentura,CPN,1/set, CS
VENT-10	FrontShroudVentura4,1/ea, CS
VENT-06	FrontTillerCoverVentura3&4,1/ea, CS
VENT-07	RearTillerCoverVentura3&4,1/ea, CS
VENT-08	ControllerVentura3&4,1/ea, CS
VENT-05	FacepltW/Switches18/20CapSeatVentura3&4, CS
VENT-02	TillerSwivelAssemblyVentura3&4,1/ea, CS
VENT-04	FacepltW/Switches18FshSeatVentura3&4/1ea, CS
VENT-03	HandlebarAssemblyVentura3&4,1/ea, CS
VENT-01	Forward/ReverseAssembly,Ventura3&4,1/ea, CS
VENT-39	ShroudDecorKitVentura4Red, CS
VENT-40	ShroudDecorKitVentura4Blue, CS